

Chapter 12

consists of two parts:

- Glossary:
explanations of machine specific terms and abbreviations.
- Key Word Index:
alphabetical listing of machine specific terms with
reference to the explanation.



Chapter 12

List of contents

1	Glossary	1
2	Key Word Index.....	9

1

ADC Compact Plus Glossary

μGy	Micro Gray – x-ray dose 1Gy= 0.87 x 10 ⁻² J/kg
A / D Converter	Analog Digital Converter ; the part that translates an analog signal into a digital signal, which can be handled by computers.
ADC	Agfa Diagnostic Center
Adjustloop	Service Program for mechanical adjustments of the CHM-module step by step cycle of a 3.5x43cm IP + Cassette, to define the correct stepper motor steps controlled by the Service PC
AE – title	Application Entity
AOS	Adonis Operating System Operating system for Agfa CPU's (Gemini, Goliath, David)
APIP	Agfa Picture Archiving Protocol
AUI	Attachment Unit Interface (External „Box“ to connect to a network)
AS	Archive Station
BOL	Begin Of Line ; Sensor for exact positioning the scan surface along the x-axis. A PIN diode is hit by the laser beam before it reaches the image plate. The Pin diode creates a digital signal, called begin of line signal. The distance (time) between this signal and the IP must be adjusted.
BOS	Begin Of Scan ;
Browser	Software that provides an interface to the World Wide Web
BSP	Boot Support Package : boot program for Agfa CPU's
Calibration	Procedure on ADC to make a homogeneous exposure look like one on film. Algorithm to counterbalance irregularities in the scanner of ADC.
CCM	Configuration and Customization Manager (Tool to edit the configuration file <i>adc.cpf</i>)
CHM	Cassette Handling Module
Collimation	Determination of the region in the image where the interesting data is at. On that region the image processing is applied..

Controller	Is connected ahead the laser recorder. The controller is responsible for the image reception from the diagnostic equipment, the image processing, layout of the images on the film, and image transmission to the machine.
cPCI	compact Peripheral Component Interconnect – bus system used in ADC Compact Plus
CPF – file	Customization Parameter File ; file that contains settings for all parameters that can be modified to configure an ADC System according to local needs.
CPU	Central Processing Unit .
CR	Computed Radiography
CU-Filter	Copper Filter, used for control of image quality
DAC	Digital Analog Converter
Decomposition	The original image is being split up in a set of bandpass-filtered images
Diagnostic Logger	Tool on ADC to keep a lot of reduced images with the corresponding image data on the HD.
DICOM	D igital I maging and CO munications in M edicine;
Dose, x-ray dose	= kV x mA x sec see also μGy = Micro Gray
Dosimeter	Device to measure x-ray dose
DR	D irect R adiography
DRA	D irect R emote A ccess (external product name: AGFATEC LINK)
Dynamic Range	Exposure range in which ADC can get usable image data.
Ethernet	A network standard for the hardware and data link levels.
Fastscan direction	
FOI	F iber O ptic I nterface. Digital interface for the transmission of image and control data via a light-leading fiber.
FSE	F ield S ervice E ngineer
FTP	F ile T ransfer P rotocol. The Internet service that transfers files from one computer to another. (Program used for transmission of files in the Internet)
Grid	to reduce scattered radiation to increase sharpness

GSC	Global Support Center
GUI	Graphical User Interface
HT power supply	High Tension power supply
HCP	Hard Copy Printer
HDD	Hard Disk Drive
Heel Effect	inhomogeneous exposure
HeNe Laser	Helium Neon Laser
HiRes	High Resolution: pixel size = 100 micrometer
HIS	Hospital Information System
Histogram	graphical display of the distribution of gray levels
HUB (Ethernet)	box to interconnect network hosts with Twisted Pair cable. Also called concentrator.
HUB (in AGFA devices)	Switchbox (Hard- or Software) which switches or multiplexes different channels, e.g. AMDI and Service Channel.
I/O BUS	Input / Output BUS. System consisting of a defined cable, cable connection, and signal for the parallel transmission of control data. The machine uses a bus system of 8 V level and 26 parallel lines.
ID Station	IDENTIFICATION Station
IMOS	IMAGE MONITOR SOFTWARE -> succession for SMA
Internet	The global computer network, composed of thousands of Wide Area Networks (WANs) and Local Area Networks (LANs), that uses TCP/IP to provide world-wide communication to homes, schools, businesses, and governments. The World Wide Web runs on the Internet.
IP	for ADC: Image Plate (phosphor plate)
IP	Image Processing
IP	(Internet Protocol) Internet software that divides data into packets for transmission over the Internet. Computers must run IP to communicate across the Internet. See also TCP.
IP-Address	(Internet Protocol Address) The standard method which identifies an internet connected computer.

Java	A general-purpose programming language created by Sun Microsystems. Java can be used to create Java applets. A Java program is downloaded from the Web server and interpreted by a program running on the machine containing the Web browser.
Java applet	A short program written in Java that is attached to a World Wide Web page and executed by the browser machine.
JavaScript	A cross-platform, World Wide Web scripting language developed by Netscape Communications. JavaScript code is inserted directly into the HTML page.
LAN	Local Area Network. Network technology, designed to connect computers over short distances. It is possible to connect the LAN with the Internet or to make a configuration into an intranet.
Laser	Light Amplification by Stimulated Emission of Radiation:
LCD	Liquid Crystal Display
LED	Light Emitting Diode; Semi conductor emitting light
Leeds Phantom	Image Plate to check the technical image quality of the ADC. The Leeds Phantom was developed for Agfa by the University of Leeds.
Low Pass Filtering	Butterworth filter in SAB – board to eliminate frequencies > 250 kHz.
LUT	Look Up Table
MD – Plate	AGFA MD → Medium Definition Image Plate
MFA	Machine Factor A Machine specific value indicating the photo multiplier sensitivity MFA is the log of the PM high voltage which results in a scan average level SAL1800 at a given dose of 1mR (~20 µGy) and a speed class 200
MFB	Machine Factor B
MI-CSO	Medical Imaging – Customer Support Operations; Agfa department in charge of world wide service.
MIMOSA	Medical Image Operating System Agfa; Agfa's workstation operating software
MODEM	MODULATOR DEMODULATOR – device to connect via telephone line to another computer.

Monitor – level	lowest software level on CPU (like BIOS on DOS systems) stored on EEPROM
mR	milli Röntgen , measure for x-ray dose 1 mR corres. 8.7 µGy
MS-Board	M ulty S upply B oard; is part of the Power Unit, detects supply voltage, controls the mechanical periphery of the digitizer and protects the stepper motor boards.
MUSICA	M ultiple S cale I mage C ontrast A mplification; Agfa's image processing software
Network location	In a URL, the unique name that identifies an Internet server. A network location has two or more parts, separated by periods, as in my.network.location. Also called host name and Internet address.
Node	Nodes in the I/O bus system = printed circuit board in the I/O bus.
Nullmodem	RS 232 (RS242) cable with crossed Transmit / Receive Line
NVE	N ame v alue f ile e ditor; editor system for parameter
OLUT	O utput L ook-Up T able
OS	O perating S ystem
Password	A text string that allows a user access to an Internet service, if the service requires it.
Photomultiplier Tube (PMT)	P hoto M ultiplier T ube: opto – electronical sensitive device to convert light (laser emission) in current and with a I / V –Converter into voltage
Pixel	Scanning point on the film or plate. A maximum of 2048 per line may be read by the Photomultiplier.
Preview Monitor	displays the image and the name of the patient. This monitor enables to roughly check whether the exposure was successful.
PRID	P REVIEW and I D Station, installed on one PC
PS	P rocessing S tation
Pyramidal Image	Image file format of images on VIPS
Quantisation	compression of image files
RAM	R andom A ccess M emory volatile main memory of computers

RAM – DISK	virtual harddisk simulated in the CPU main memory Gemini RAM-Disk contains all the machine dependent parameters, e.g. stepper motor steps, info counters
REM Tool	Debugging tool for AOS logfiles
Reset	Machine reset into a defined machine status. Various checking routines are carried out during a reset.
RIS	R adiology I nformation S ystem
ROI	R egion O f Interest
ROM	R ead O nly M emory
RS232 interface	Serial interface which converts computer internal parallel information into serial bits, and vice versa.
SAL	S can A verage L evel Digitized photo multiplier signal of an average of several hundred scanned lines Range SAL0 – SAL4095 (12bit)
Scan Position Adjustment	
SCP	S OFT C OPY Tool
SCSI	S mall C omputer S ystems I nterface – interface to connect peripherals to computers (e.g. HDD, CD-ROM)
Server	A computer that shares its resources, such as printers and files, with other computers on the network.
Service Host	
Service interface	RS232 standard interface for the connection of the Service PC.
Shading Calibration	position dependent sensitivity calibration; calibration of each pixel in a line
Shell	User Interface of the VME – AOS
Slowsan direction	Transport Direction of Image Plate through the scan unit. The stepper motor speed is calculated such that the distance between two scan lines equals the distance between two pixels in a line.
Speed Class (SC)	Dose=1mR + SC100 => D=1 above fog + base A film-screen-system with speed class 100 which is exposed with a dose of 1mR results in a density 1 above fog+base.

Square Root Compression	Signal $\sim \sqrt{N}$. Method to quantize Signal in 12 similar block sizes
Standard Res	Standard Resolution: pixel size = 150 micrometer
Subnet number	A part of the internet address which designates a subnet
TCP	T ransmission C ontrol P rotocol. Internet networking software that controls the transmission of packets of data over the Internet. Among its tasks, TCP checks for lost packets, puts the data from multiple packets into the correct order, and requests that missing or damaged packets be resent. Computers must run TCP to communicate with World Wide Web servers.
TCP/IP	T ransport C ontrol P rotocol / I nternet P rotocol. Common Protocol for Networks. Used in the Internet
Testsheet	Image plate with testpattern, print to check the geometric quality of ADC:
Text Field	Part of the film displaying patient, hospital and image processing information
UI	U ser I nterface; part of a computer program that handles interaction with the user.
UPS	U ninterruptible P ower S upply
URL	U niform R esource L ocator; address code for multimedia documents.
UTP - cable	U nshielded T wisted P air (Network transfer medium) cable. Cable to connect to a computer network.
VIPS	V iewing and I mage P rocessing S tation
VME bus	V ersa M odule E urope b us Bus system from Motorola (used in ADC Digitizers) Bus for the parallel image data transmission. 16 bit bus.
W/L	W indow and L eveling of digital images.
X-Rite densitometer	Automatic density measuring instrument manufactured by X-Rite.





2 Key Word Index

A

Access to Service Level 3.2/3

Accessories 11/4

Accessory Kit 1/6

ADC Cassettes

ADC MD30 Image Plate Principle 2/6

Asymmetric Cassette 2/7

Mammography Cassette 2/6

Standard Cassette 2/5

Adjustments

Overview 3.6/1

Adjustments, mechanical

Adjustment Cycle 3.6/8

Guiding Plate in Postscan Unit 3.6/4

Rotation Unit 3.6/2

Suction Cups 3.6/7

Adjustments, optical

Shading Calibration 3.6/11

B

Backup 1/21

Banding 1/25

BOL (Begin Of Line) 2/47

BOS (Begin Of Scan) 2/47

C

Calibration 1/25

Calibration Principle

BOL (Begin Of Line) 2/47

BOL/BOS Adjustment 2/49

BOS (Begin Of Scan) 2/47

Definition 2/47

Scan Position Adjustment 2/49

Shading Calibration 2/49

Cassette

Putting out 2/40

Transporting to Output 2/40

CCM Tool 1/16, 3.2/1

Certificates 1/3

Checking the Image Quality

Evaluation of a Flatfield 9/12

Exposure of a Flatfield 9/10

Checklist

Installation 1/26

Installation Planning 11/12

Codings

GS IOB-SIN-5Step Boards 3.4/1

GS IOB-SIN-Step Boards 3.4/1

Conditions, ambient

Emissions 11/8

Floor Conditions 11/8

Humidity 11/8

Light Tightness 11/8

Magnetic Fields 11/8

Temperature 11/8

Connections, electrical 11/9

Connection Cables 11/9

Fuse Protection 11/9

Network Connection 11/9

Power Connection 11/9

Connectors 1/14

Controls 1/14

cPCI Power Supply 2/19

CPU-Ariel 2/16

Cycle of a Cassette in the Digitizer

Cassette Unit 2/39

Cassette Unit 1/2 2/31

Input Buffer 2/30

IP Transport Postscan/Erasure Unit 2/36

IP Transport Prescan 2/32

Output Buffer 2/40

Rotation Unit 2/34

Scan Unit 2/34

D

Diagrams 4/27

Digitizer

- ADC System Components 2/1
- Work Flow of the ADC System 2/4

Dimensions 11/6, 11/11

E

Electrical Data

- Frequency 11/10
- Leakage current 11/10
- Power consumption 11/10
- Rated Voltage 11/10

Erase Control Board 2/12

F

Flatfield 1/22, 3.6/16

Form for Network Parameters 11/14

Fuses 3.4/2

H

Hard disk 2/20

I

Identification Diagrams 4/1

Installation

- Checklist 1/26
- Optic Module 3.5/13
- Planning Checklist 11/12
- Preparations before 1/5
- Requirements 1/16
- Start-up Procedure 1/19

Interlock switch 3.1/4

IP

- Positioning 2/34
- Putting back into the Cassette 2/38
- Scanning 2/35
- Taking out 2/33
- Transporting to the Erasure Unit 2/37
- Transporting to the Scanner 2/33
- Unloading from Scan Unit 2/35

L

Laser Diodes 2/26

M

Maintenance

- Checking Electrical Functions 9/8
- Recommended Spare Parts 9/3
- Required Spare Parts 9/3
- Required Tools/Material 9/2

Maintenance Checklist

- ADC System 9/14
- Digitizer 9/13

Maintenance Points

- Cassette Unit 9/4
- cPCI-Rack 9/6, 9/8
- Erasure Unit 9/7
- Erasure Unit Fan 9/5
- Image Plates/Cassettes 9/9
- Power Unit 9/4, 9/8
- Scan Unit 9/6
- Transport Units 9/5

Modular Structure of the Digitizer

- Cassette Unit 1/2 2/22
- cPCI-Rack 2/14
- Erasure Unit 2/28
- Input Buffer 2/9
- IP Transport Unit Postscan 2/29
- IP Transport Unit Prescan 2/23
- Output Buffer 2/13
- Overview of Modules 2/8
- Power Unit 2/10
- Rotation Unit 2/24
- Scan Unit 2/25

Multi Supply Board 2/11

O

Operation Terminal 1/15

Optic Module Housing 2/26

Overview of ADC Web Pages 3.2/9

P

Photomultiplier (PMT) 2/26

Polygon 2/26

**R****Recommended Tools 3.2/1****Remote Service 3.2/7****Removing**

- Optic Module 3.5/10
- Scan Unit 3.5/14
- Shipping Brackets 1/7

Replacement

- 5fold Stepper Motor Board 3.5/9
- Cassette Unit 3.5/15
- Erasure Unit 3.5/1
- Erasure Unit Fan 3.5/3
- Lamps 3.5/2
- Optic Module 3.5/10
- Photomultiplier Tube (PMT) 3.5/8
- Scan Unit 3.5/8
- Stepper Motor Boards of the Robot 3.5/7
- Transport Units 3.5/4
- Transport Units (Pre- and Postscan) 3.5/6
- Vacuum Pump 3.5/4

S**Safety Instructions 1/1, 3.1/1, 11/1**

- Cassette Unit 1/2, 3.1/2, 11/2
- Cassettes 3.1/3
- Laser Products 1/2, 3.1/2, 11/2

Safety Notes 3.1/3, 3.5/1**Safety Regulations**

- Certificates 11/3

Safety Switch 3.1/4**SAL (Scan Average Level) 2/49****Scan Master Board 2/17****Scope of Delivery 11/4****Service Key 3.1/4****Service PC 3.2/2****Service Program**

Overview 3.2/5

Shipping Brackets 1/7**Signal and Data Flow**

- Data Flow 2/47
- Signal Flow 2/45

Software Architecture

- Diagnostics Software 2/42
- Hardware Modification Floppy 2/43
- Web Pages 2/43

System Check 1/21**System Overview 11/5****T****Test Equipment 3.2/1****Test Sheet 1/23****Transport Path 11/6****Transport Robot 2/32****Troubleshooting**

- Cassette Unit 3.3/3
- Erasure Unit 3.3/12
- Input Buffer 3.3/2
- Overview 3.3/1
- Postscan Transport Unit 3.3/14
- Prescan Transport Unit 3.3/6
- Rotation Unit 3.3/5
- Scanner 3.3/11

V**Voltage Adjustment 1/18****W****Warning Labels**

- Erasure Unit 3.1/7
- Input Buffer 3.1/6

Weights 11/11

